


# PLACA BLOQUEADA LPS T TÍBIA MEDIAL EM TITÂNIO

*LOCKING PLATE LPS  
T MEDIAL TÍBIA IN TITANIUM*

*PLACA BLOQUEADA LPS T TIBIA  
MEDIAL EN TITANIO*



ORTOBIO®



# PLACA BLOQUEADA LPS T TÍBIA MEDIAL EM TITÂNIO

LOCKING PLATE LPS T MEDIAL TIBIA IN TITANIUM

PLACA BLOQUEADA LPS T TIBIA MEDIAL EN TITANIO




Tem como indicação, finalidade e/ou uso a estabilização e fixação de diversos tipos de fraturas, proporcionando a recuperação de fraturas na Tíbia Medial. A conexão rígida entre os parafusos de bloqueio e a placa permite melhores resultados clínicos em fraturas no platô tibial, fraturas simples e cominutas e em ossos osteoporóticos da região medial da tíbia.



Its indication, purpose and/or use is the stabilization and fixation of different types of fractures, providing the recovery in the medial tibia. The rigid connection between the locking screws and the plate allows better clinical results in simple and comminuted fractures, and in osteoporotic bones in the medial region of the tibia.



Su indicación, finalidad y/o uso es la estabilización y fijación de diferentes tipos de fracturas, proporcionando la recuperación de las fracturas en la tibia medial. La conexión rígida entre los tornillos de bloqueo y la placa permite mejores resultados clínicos en las fracturas simples y conminutas y en huesos osteoporóticos en la región medial de la tibia.



ANVISA 80062900031

**Titânio Puro** (NBR ISO 5832-2)  
**Pure Titanium** (NBR ISO 5832-2)  
**Titanio Puro** (NBR ISO 5832-2)



**Biologicamente compatível,**  
evitando reação inflamatória.

**Biocompatible,** avoiding  
inflammatory reactions.

**Biologicamente compatible,**  
evitando reacción inflamatoria.



**Formato anatômico**  
**Anatomic shape**

**Formato anatómico**



**Acabamento polido**  
**Polished finish**

**Acabado pulido**



**Furos de bloqueio:**  $\varnothing 5,5\text{mm}$ ,  
**e simples:**  $\varnothing 5,5\text{mm}$

**Locking holes:**  $\varnothing 5,5\text{mm}$ ,  
**and simple:**  $\varnothing 5,5\text{mm}$

**Perforaciones de bloqueo:**  
 $\varnothing 5,5\text{mm}$  **y simples:**  $\varnothing 5,5\text{mm}$



**Utilizado com parafusos:**

Bloqueio:  $\varnothing 5,0\text{mm}$   
Cortical:  $\varnothing 4,5\text{mm}$   
Esponjoso:  $\varnothing 6,5\text{mm}$

**Use with screws:**

Locking:  $\varnothing 5,0\text{mm}$   
Cortical:  $\varnothing 4,5\text{mm}$   
Cancellous:  $\varnothing 6,5\text{mm}$


**Utilice con tornillos:**

De bloqueo:  $\varnothing 5,0\text{mm}$   
Cortical:  $\varnothing 4,5\text{mm}$   
Esponjoso:  $\varnothing 6,5\text{mm}$

# PLACA BLOQUEADA LPS T TÍBIA MEDIAL EM TITÂNIO

LOCKING PLATE LPS T MEDIAL TIBIA  
IN TITANIUM

PLACA BLOQUEADA LPS T TIBIA  
MEDIAL EN TITANIO



ANVISA 80062900031

<b>CÓDIGO</b> <i>CODE</i> <b>CÓDIGO</b>	<b>Nº FUROS</b> <i>Nº HOLES</i> <b>Nº PERF.</b>	<b>LARGURA (MM)</b> <i>WIDTH (MM)</i> <b>ANCHO (MM)</b>	<b>COMPRIM. (MM)</b> <i>LENGHT (MM)</i> <b>LONGITUD (MM)</b>	<b>ESPESSURA (MM)</b> <i>THICKNESS (MM)</i> <b>ESPESOR (MM)</b>
5713-103-002	2	14	64	3,0
5713-103-003	3	14	84	3,0
5713-103-004	4	14	104	3,0
5713-103-005	5	14	124	3,0
5713-103-006	6	14	144	3,0
5713-103-007	7	14	164	3,0
5713-103-008	8	14	184	3,0
5713-103-009	9	14	204	3,0
5713-103-010	10	14	224	3,0
5713-103-011	11	14	244	3,0
5713-103-012	12	14	264	3,0